Points to Practice

DISINFECT

- Disinfect the silicone septum with an appropriate antiseptic agent.
- Disinfect prior to each access.

ACCESS

Syringe:

- Access the BD Q-Syte device with an ISO luer lock or luer slip syringe.
- Connect using a straight-on approach.

Administration Set:

- When using an administration set with a spin nut, insert the slip luer into the BD Q-Syte device. Then secure the connection with the spin nut.

Note: Do not use needles or blunt cannula with the BD Q-Syte device.



Blood Withdrawal:

- Use ISO syringe or BD Luer-Lok™ access device.

FLUSH

- Flush according to facility policy or until device is clear. If unable to clear replacement is recommended.
- To minimize reflux of blood into the vascular access device, it is recommended to clamp before disconnecting.
- When disconnecting, hold the BD Q-Syte device firmly in place to prevent accidental removal.

The BD Q-Syte device has a unique luer penetrating septum design that provides higher flow rates, fluid path visualization, and consistent dependable performance.

Non-tortuous and unrestricted fluid path

Helping all people
live healthy lives

NOTE: Refer to package insert for precautions and complete instructions for use.

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
⁵ Data on file
⁶ <http://www.cms.hhs.gov/apps/media/press/release.asp>
⁷ Centers for Medicare & Medicaid Services. (2008, April 14) CMS proposes additions to list of hospital-acquired conditions for fiscal year 2009. Available at: http://www.cms.hhs.gov/apps/media/fact_sheets.asp

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¹ Rupp ME, Sholtz LA, Jourdan DR, et al. Outbreak of bloodstream infection temporally associated with the use of an intravascular needleless valve. *CID*. 2007;44:1408-1414.
² Salgado CD, Chinnes L, Paczesny TH, Cantej JR. Increased rate of catheter-related bloodstream infection associated with use of a needleless mechanical valve device at a long-term acute care hospital. *Infect Control Hosp Epidemiol*. 2007;28:684-688.
³ Klevens RM, Edwards JR, Richards CL, et al. Estimating health care-associated infections and deaths in U. S. hospitals, 2002. *Public Health Reports*. 2007; 122:160-166.
⁴ Karchmer TB, Wood C, Ohl CA, et al. Contamination of mechanical valve needleless devices may contribute to catheter-related bloodstream infections. *SHEA 2006 Presentation Number: 221 Poster Board Number: 47*.